

POSTURE CORREX



The XBACK POSTURE CORREX is an LSO brace that provides coverage from L1 - L5 incorporating a unique easy to adjust strapping system in order to correct postural deficits and acute or longstanding disorders of the lumbar, thoracic and cervical spine as well as shoulder girdle.

MODEL T106 L0627 PDAC APPROVED



XBACK FITTING & CARE INSTRUCTIONS

Simply Innovative Braces™

POSTURE CORREX MODEL T106

FOLLOW THESE INSTRUCTIONS FOR WEARING AND CARING FOR YOUR XBACK BRACE

Steps for putting on and adjusting your brace:



Remove brace from package and undo straps.



Center brace on the back with straps at top of back panel



Close the two halves of the front low and snug on the abdomen.



Take the left strap with your left hand



and toss it over the left shoulder



Reach behind you with the right hand and grasp the loose end of the strap



Pull the left strap across your back



and under your right arm / shoulder



Attach the strap to the front panel Repeat the same steps in reverse with the right strap



Now loosen both of the straps, stand up straight rotate your shoulders and shoulder blades back and down.



Pull on the straps simultaneously to the desired level of support and re-anchor the straps to the front panel. Finally you can adjust the shoulder pads for comfort

Care Instructions:

Hand wash only in cool water with a mild liquid soap and lay flat to air dry. When drying do not expose to direct heat or sunlight. Do not dry clean. If your brace has a rigid shell, the shell may be cleaned with a mild topical cleaner. Do not use strong detergents or bleaches as both may damage the brace materials.

Warranty:

All XBACK braces are warranted to be free of material or workmanship defects for 5 years. Any defective product will be repaired or replaced at the discretion of XBack Bracing Services, Inc. after review of the defect by Xback Bracing Services, Inc. There are no other implied or expressed warranties.



Step 1: Strap is brought under the arm

Step 2: Strap is then brought over the shoulder

Step 3: strap is brought across the back

Step 4: Strap is brought under the opposite arm/shoulder

Step 5: Strap is attached to the front panel

Step 6: Repeat in the opposite direction for other strap.



XBack Bracing Services, Inc. P.O. Box 100 / 341A West Main Street, Birdsboro, PA 19508, USA
610-404-4900 - Fax 610-404-4905 www.XBackBrace.com / info@XBackBrace.com

